

BEAUTY.GIENE

HOSPITAL GRADE DISINFECTANT

HY.GIENE
AUSTRALIA

*Multipurpose cleaning in
hair, beauty & lifestyle
industries*

PRODUCT FEATURES & BENEFITS

- Multipurpose cleaner, degreaser, sanitiser and hospital grade disinfectant that kills 99.99% of germs.

PRODUCT USE

Mix BEAUTY.GIENE with water at

DILUTION RATE 1:4 (20%)

100ml of BEAUTY.GIENE into 500ml spray bottle.

Spray & Wipe: Spray onto surface and allow to stand for 10 minutes. Wipe using a clean damp cloth and allow to dry.

Only suitable as a hard surface disinfectant - stainless steel, bathrooms, chairs and equipment.

Mopping: use BEAUTY.GIENE at a dilution of 1:80 (1/4 cup in a bucket). if oiling is very heavy and disinfection is critical, clean first with BEAUTY.GIENE, then use a fresh solution afterwards to disinfect.

Soaking Solution: 1:5 (200ml water) for combs, clippers & scissors.

Do not mix with soaps, detergents or other chemicals or deactivation may occur.

PRODUCT AVAILABILITY

1L product code 80101

5L product code 80105

15L product code 80120



TECHNICAL DATA

Colour: Thin, clear, red-coloured liquid

Odour: Characteristic (unscented)

pH: 11.0 to 11.5

SAFETY INFORMATION

Refer to product label and Safety Data Sheet before using this product.

This sheet contains information that you may need to protect your employees and customers against health or safety hazards associated with our product.

Not a Dangerous Good according to ADG Code, IATA or IMDG/IMSBC criteria.

ISSUE DATE: FEBRUARY 2025